

## **מעבדת תרומת דם**

מגישות:

אורית רוקח 213279631

שירהל מיתר 214439101

## תוכן העניינים

3..... מערכת: תיאור מילולי

### שלב א

4..... ERD של המערכת

5..... DSD של המערכת

6..... CreateTable

10..... DropTable

11..... Insert

12..... SelectAll

13..... Desc

15..... DataGenerators

17..... גיבוי

### שלב ב

19..... שאלות בחירה

23..... שאלות עידכון

25..... שאלות מחיקה

27..... שאלות עם פרמטרים

31..... אילוצים

### שלב ג

32..... Procedures

36..... Functions


38..... Mains

## **תיאור מילולי:**


**שם המערכת:** מעבדת תרומות דם.

**תיאור המערכת:** המערכת היא מערכת לניהול מידע על מכשירים, בדיקות רפואיות, תרומות דם וציווד רפואי. הנתונים המרכזיים במערכת כוללים מידע על מכשירים רפואיים (כמו מספר סידורי, סוג, יצרן ותאריכי רכישה ותחזוקה אחרונה), בדיקות רפואיות (כולל סוג בדיקה, תאריך ותוצאות), תרומות דם (כולל סוג דם, תאריך התרומה וכמות דם) וציווד רפואי נלווה (כולל תאריכי רכישה ותחזוקה, דגם ומיקום). המערכת מטרתה לאפשר ניהול יעיל ומדויק של כלל המידע הרלוונטי, כולל הקשרים בין הנתונים, לצורך תמיכה בתהליכי רפואיים, ניהול מלאי וציווד, וניטור וניהול תרומות דם.

## ERD של המערכת:



## DSD של המערכת



## CreateTable

```
CREATE TABLE Device1
(
    DeviceId varchar2(50) NOT NULL,
    DeviceType INT NOT NULL,
    SeriesNum INT NOT NULL,
    VersionD INT NOT NULL,
    Manufacturer varchar2(50),
    PurchaseDate DATE,
    LastMaintenanceDate DATE,
    LocationD varchar2(50),
    PRIMARY KEY (DeviceId)
);
```

```
CREATE TABLE Test_Tube1—
(
    TestTubeID varchar2(50) NOT NULL,
    Material INT,
    VolumeCapacity INT,
    PurchaseDate DATE,
    ExpiryDate DATE,
    LocationTestTube varchar2(50),
    PRIMARY KEY (TestTubeID)
);
```

```
CREATE TABLE Test_1—
(
    TestID varchar2(50) NOT NULL,
```

```

MainTestType INT NOT NULL,
SubTestType INT NOT NULL,
TestDate DATE NOT NULL,
ResultT varchar2(50) NOT NULL,
BloodDonationID varchar2(50) NOT NULL,
PRIMARY KEY (TestID),
UNIQUE (TestDate),
UNIQUE (BloodDonationID)
);

```

```

CREATE TABLE Test_Device1—
(
TestID varchar2(50) NOT NULL,
DeviceId varchar2(50) NOT NULL,
PRIMARY KEY (TestID, DeviceId),
FOREIGN KEY (TestID) REFERENCES Test_1(TestID),
FOREIGN KEY (DeviceId) REFERENCES Device1(DeviceId)) ;

```

```

CREATE TABLE DeviceTest_Tube1—
(
DeviceId varchar2(50) NOT NULL,
TestTubeID varchar2(50) NOT NULL,
PRIMARY KEY (DeviceId, TestTubeID),
FOREIGN KEY (DeviceId) REFERENCES Device1(DeviceId),
FOREIGN KEY (TestTubeID) REFERENCES Test_Tube1(TestTubeID)
);

```

```

CREATE TABLE Blood_Donation1

```

```
(
    DonationID varchar2(50) NOT NULL,
    DonorID varchar2(50) NOT NULL,
    BloodType varchar2(50) NOT NULL,
    DonationDate DATE,
    VolumeDonated INT,
    TestResultID varchar2(50),
    TestID varchar2(50) NOT NULL,
    PRIMARY KEY (DonationID),
    FOREIGN KEY (TestID) REFERENCES Test_1(TestID),
    UNIQUE (DonorID)
);
```

CREATE TABLE Stove1—

```
(
    StoveId varchar2(50) NOT NULL,
    SeriesNum INT NOT NULL,
    VersionStove INT NOT NULL,
    Manufacturer INT,
    PurchaseDate DATE NOT NULL,
    LastMaintenanceDate DATE NOT NULL,
    LocationStove varchar2(50) NOT NULL,
    DeviceId varchar2(50) NOT NULL,
    PRIMARY KEY (StoveId),
    FOREIGN KEY (DeviceId) REFERENCES Device1(DeviceId)
);
```

CREATE TABLE Refrigerator1—



```
(  
    RefrigeratorId varchar2(50) NOT NULL,  
    ModelR INT,  
    PurchaseDate DATE,  
    LastMaintenanceDate DATE,  
    LocationR varchar2(50),  
    TestTubeID varchar2(50) NOT NULL,  
    DonationID varchar2(50) NOT NULL,  
    PRIMARY KEY (RefrigeratorId),  
    FOREIGN KEY (TestTubeID) REFERENCES Test_Tube1(TestTubeID),  
    FOREIGN KEY (DonationID) REFERENCES Blood_Donation1(DonationID)  
);
```

## **DropTable**

```
;DROP TABLE Refrigerator1  
;DROP TABLE Stove1  
;DROP TABLE Blood_Donation1  
;DROP TABLE DeviceTest_Tube1  
;DROP TABLE Test_Device1  
;DROP TABLE Test_1  
;DROP TABLE Test_Tube1  
;DROP TABLE Device1
```

## Insert

זה כמות אדירה ולכן נרשום רק כמה שורות...

```
INSERT INTO Device1 (DeviceId, DeviceType, SeriesNum, VersionD,  
Manufacturer, PurchaseDate, LastMaintenanceDate, LocationD) VALUES  
('YvSGkO6l6fRwXlaf5But', 31, 760, 9, 'PvsL2vp1kgjRPSYTZaTE', TO_DATE('2003-  
09-16', 'YYYY-MM-DD'), TO_DATE('2010-04-13', 'YYYY-MM-DD'),  
'evqI5CvnJlWT5sjpFOr8');
```

```
INSERT INTO Device1 (DeviceId, DeviceType, SeriesNum, VersionD,  
Manufacturer, PurchaseDate, LastMaintenanceDate, LocationD) VALUES  
('eR4LgsUjk1PKOmGQkmDp', 46, 436, 2, 'WiFxVJLbwrlcYzZdpL5G',  
TO_DATE('2010-04-04', 'YYYY-MM-DD'), TO_DATE('2010-09-25', 'YYYY-MM-DD'),  
'mFjXIEDHUPMHZDlvApsT');
```

```
INSERT INTO Device1 (DeviceId, DeviceType, SeriesNum, VersionD,  
Manufacturer, PurchaseDate, LastMaintenanceDate, LocationD) VALUES  
('3h5ggjXUeEB2Gr8WzzWa', 73, 748, 9, 'sApvxdDVio2a9G4AMJ3C',  
TO_DATE('2020-03-12', 'YYYY-MM-DD'), TO_DATE('2000-04-20', 'YYYY-MM-DD'),  
'r964gYNL8H8739jrYgrf');
```

### SelectAll

SELECT\* from Device1;

SELECT\* from Test\_Tube1;

SELECT\* from Test\_1;

SELECT\* from Test\_Device1;

SELECT\* from DeviceTest\_Tube1;

SELECT\* from Blood\_Donation1;

SELECT\* from Stove1;

SELECT\* from Refrigerator1;

## Desc

```
VOLUMEDONATED INTEGER Y  
TESTRESULTID VARCHAR2 (50) Y  
TESTID VARCHAR2 (50)
```

SQL> desc Stove1

| Name                | Type          | Nullable | Default | Comments |
|---------------------|---------------|----------|---------|----------|
| STOVEID             | VARCHAR2 (50) |          |         |          |
| SERIESNUM           | INTEGER       |          |         |          |
| VERSIONSTOVE        | INTEGER       |          |         |          |
| MANUFACTURER        | INTEGER       | Y        |         |          |
| PURCHASEDATE        | DATE          |          |         |          |
| LASTMAINTENANCEDATE | DATE          |          |         |          |
| LOCATIONSTOVE       | VARCHAR2 (50) |          |         |          |
| DEVICEID            | VARCHAR2 (50) |          |         |          |

SQL> desc Refrigerator1

| Name                | Type          | Nullable | Default | Comments |
|---------------------|---------------|----------|---------|----------|
| REFRIGERATORID      | VARCHAR2 (50) |          |         |          |
| MODELNR             | INTEGER       | Y        |         |          |
| PURCHASEDATE        | DATE          | Y        |         |          |
| LASTMAINTENANCEDATE | DATE          | Y        |         |          |
| LOCATIONNR          | VARCHAR2 (50) | Y        |         |          |
| TESTTUBEID          | VARCHAR2 (50) |          |         |          |
| DONATIONID          | VARCHAR2 (50) |          |         |          |

SQL> |

SQL> desc Test\_Device1

| Name     | Type          | Nullable | Default | Comments |
|----------|---------------|----------|---------|----------|
| TESTID   | VARCHAR2 (50) |          |         |          |
| DEVICEID | VARCHAR2 (50) |          |         |          |

SQL> desc DeviceTest\_Tube1

| Name       | Type          | Nullable | Default | Comments |
|------------|---------------|----------|---------|----------|
| DEVICEID   | VARCHAR2 (50) |          |         |          |
| TESTTUBEID | VARCHAR2 (50) |          |         |          |

SQL> desc Blood\_Donation1

| Name          | Type          | Nullable | Default | Comments |
|---------------|---------------|----------|---------|----------|
| DONATIONID    | VARCHAR2 (50) |          |         |          |
| DONORID       | VARCHAR2 (50) |          |         |          |
| BLOODTYPE     | VARCHAR2 (50) |          |         |          |
| DONATIONDATE  | DATE          | Y        |         |          |
| VOLUMEDONATED | INTEGER       | Y        |         |          |
| TESTRESULTID  | VARCHAR2 (50) | Y        |         |          |
| TESTID        | VARCHAR2 (50) |          |         |          |

SQL> desc Stove1

| Name         | Type          | Nullable | Default | Comments |
|--------------|---------------|----------|---------|----------|
| STOVEID      | VARCHAR2 (50) |          |         |          |
| SERIESNUM    | INTEGER       |          |         |          |
| VERSIONSTOVE | INTEGER       |          |         |          |
| MANUFACTURER | INTEGER       | Y        |         |          |

```

SQL> desc device1
Name          Type          Nullable Default Comments
-----
DEVICEID      VARCHAR2 (50)
DEVICETYPE    INTEGER
SERIESNUM     INTEGER
VERSIOND      INTEGER
MANUFACTURER  VARCHAR2 (50) Y
PURCHASEDATE  DATE          Y
LASTMAINTENANCEDATE DATE        Y
LOCATIOND       VARCHAR2 (50) Y


SQL> desc Test_Tube1
Name          Type          Nullable Default Comments
-----
TESTTUBEID    VARCHAR2 (50)
MATERIAL      INTEGER      Y
VOLUMECAPACITY INTEGER      Y
PURCHASEDATE  DATE          Y
EXPIRYDATE    DATE          Y
LOCATIONTESTTUBE VARCHAR2 (50) Y

SQL> desc Test_1
Name          Type          Nullable Default Comments
-----
TESTID        VARCHAR2 (50)
MAINTESTTYPE  INTEGER
SUBTESTTYPE   INTEGER
TESTDATE      DATE

```

# Data Generators









\*על ידי תוכנית בC++



\*על ידי קצבי exel



|    | A            | B         | C           | D           | E          | F           | G            | H            |
|----|--------------|-----------|-------------|-------------|------------|-------------|--------------|--------------|
| 1  | StoveId      | SeriesNum | VersionStov | Manufacture | PurchaseDa | LastMainter | LocationSto  | DeviceId     |
| 2  | OJ8Wfwytgs   | 174       | 6           | u61GURC3Q   | 18/05/2012 | 13/08/2011  | BoPRzd0Wq    | 3JwlU9dGEK   |
| 3  | PhjEe2RBjeZ  | 558       | 4           | DrP9PLdrDE  | 11/03/2000 | 26/02/2017  | VutQ1URAbg   | GJuUmCFHC    |
| 4  | ucv1cbpZpG   | 450       | 6           | ixWGvNwzjV  | 09/08/2021 | 30/05/2008  | JLVnsJ16iwT  | VeaALT3zZXY  |
| 5  | OtW1RvJk2c   | 232       | 1           | csnDyAO7cT  | 20/05/2014 | 18/03/2022  | OGHZcYN8ey   | c2S8bsiiL3b  |
| 6  | qmbzvsh75    | 619       | 6           | bNMtNT0Hb   | 18/01/2009 | 30/11/2011  | dVLRZGP62l   | aG8voTCFXY   |
| 7  | Eh2froCbIfP  | 467       | 8           | AdEhesNB2f  | 15/12/2008 | 17/12/2018  | egFJsZF4KM   | oXmV4VXME    |
| 8  | geg5hj1XtJQ  | 785       | 2           | ohntDbjHMc  | 04/01/2015 | 03/11/2002  | yvs03aKBXv   | pOmuSefulv   |
| 9  | cAmrtwNx9l   | 348       | 10          | QtMFaAOaQ   | 20/10/2016 | 02/06/2004  | 43a7gSo2YO   | 6TZmHHpk2    |
| 10 | xP7xvGst9Jz  | 964       | 8           | xPEVVB4M6   | 14/02/2017 | 30/06/2021  | m8akuXm3u    | TDrY6r5Ed1S  |
| 11 | PqhxYb34y    | 70        | 7           | 7rcgEJWFDX  | 18/10/2021 | 22/07/2005  | dNOxhXQWJ    | VvQmJJqAG    |
| 12 | loWYjCPkxzc  | 599       | 4           | xPEVVB4M6   | 07/07/2017 | 06/07/2001  | kPJ6wn5HHf   | iEvYsrJUPTxe |
| 13 | EB9ezF823d   | 79        | 2           | mXOT9piWC   | 11/12/2009 | 08/11/2022  | FxlaHG0114   | 6TZmHHpk2    |
| 14 | fep9TtSfwdj  | 660       | 4           | x2PuvzbtSyC | 24/02/2002 | 06/04/2005  | 59dlwta97K   | VvQmJJqAG    |
| 15 | 6ro8MKCON    | 309       | 3           | RoTCaXQL8E  | 08/07/2019 | 07/07/2022  | dcCBE3VjI7n  | zVz4eeNfme   |
| 16 | UKHpG0XGv    | 604       | 7           | CZtZyasZln6 | 23/09/2011 | 04/07/2017  | n4qlODAzP6   | Dt7CfR1315C  |
| 17 | UdOKkkAQeI   | 99        | 2           | u61GURC3Q   | 02/08/2021 | 31/10/2015  | IreuSnbbyNk  | Q9AlgKi2VQI  |
| 18 | gecyR2hzohc  | 371       | 7           | h6Qa6DmKc   | 06/06/2001 | 03/06/2017  | UJdUtdHC16   | gcZRRWdAL8   |
| 19 | 6l7Ov4qMr7   | 437       | 7           | zZ7AzXqFit3 | 18/08/2023 | 19/02/2003  | rtJWWFXq3    | UpHsnKtDiV   |
| 20 | FbKBuDvdM    | 310       | 8           | RVi0pauvdW  | 18/06/2016 | 09/06/2005  | zraIB27hQn   | MHRvIrODdE   |
| 21 | pLBOWESXQ    | 310       | 2           | CZtZyasZln6 | 23/07/2020 | 29/03/2015  | Zot8YWkgxc   | YFvzR6EXWS   |
| 22 | ten6PvTviRV  | 304       | 6           | NxpUebcl3gC | 02/12/2012 | 18/08/2007  | BcwsWb2Lw    | eneDROYPYF   |
| 23 | YCPYnvqT1F   | 894       | 9           | rbCuTVjFSU5 | 07/01/2005 | 21/07/2013  | JclUslbia4og | X7PapIqOITh  |
| 24 | lBOHfxfnSfgM | 59        | 3           | Y8f1XoK7d8  | 04/03/2014 | 13/12/2023  | xG4dvyyiqR5  | cD7b4H2mT    |
| 25 | HPojJsJF6tVl | 433       | 6           | NxpUebcl3gC | 10/07/2005 | 22/02/2018  | sRXlLJKockR  | BYsYqbbqkNlJ |

## \*על ידי Data Generaor

| Name  | Date modified    | Type    | Size |
|---|------------------|---------|------|
|  DataGeneratorBloodDonation.gd     | 25/05/2024 21:48 | GD File | 1 KB |
|  DataGeneratorDevice_Test_Tube1.gd | 25/05/2024 21:37 | GD File | 1 KB |
|  DataGeneratorDevice1.gd           | 25/05/2024 21:37 | GD File | 1 KB |
|  DataGeneratorRefrigerator1.gd     | 25/05/2024 21:48 | GD File | 1 KB |
|  DataGeneratorStove1.gd            | 25/05/2024 21:44 | GD File | 1 KB |
|  DataGeneratorTest_Device1.gd      | 25/05/2024 21:48 | GD File | 1 KB |
|  dataGeneratorTest_Tube1.gd        | 25/05/2024 21:37 | GD File | 1 KB |
|  DataGeneratorTest1.gd             | 25/05/2024 21:48 | GD File | 1 KB |



## גיבוי



Source Objects Data Other

Connections CreateTable.sql DropTables.sql InsertTable.sql SelectAll.sql desc.pdc Import Tables Edit data STOVE1@LABDBWIN PL/SQL Clipboard

SQL Output Statistics

```
select t.*, t.rowid from STOVE1 t
```

Objects Files Templates Window List

<CURRENT USER>

All objects

Enter filter text...

DBMS\_Jobs

Queues

Queue tables

Libraries

Directories

Tables

BLOOD\_DONATION1

DEVICE1

DEVICETEST\_TUBE1

REFRIGERATOR1

STOVE1

TEST\_1

TEST\_DEVICE1

TEST\_TUBE1

Indexes

Constraints

Views

| STOVEID | SERIESNUM   | VERSIONSTOVE | MANUFACTURER | PURCHASEDATE             | LASTMAINTENANCE |
|---------|---|--------------|--------------|--------------------------|-----------------|
| 1       | 1bQv9CpTWWFL9wK7NvI                               | 25           | 8            | 08/07/2019               | 27/04/2013      |
| 2       | dkq829bM475ih275oNh09nQb149hKT716Mo63N1579aS      | 771          | 29           | 144 13/10/2208 17:50:21  | 10/11/2028 08:  |
| 3       | wMD693dkV856bHf72aAe575R25vMC781tT500eR0681fN     | 795          | 50           | 1459 13/01/2524 15:54:17 | 11/03/2303 16:  |
| 4       | YD773cAw787eOy719hG868gWB897bPA976kQs571nOG55yL   | 1170         | 97           | 578 11/02/2100 19:43:23  | 17/10/2396 12:  |
| 5       | nKL628Cn816aGP799bRL287gY0270yXk669kZV758aJb989pO | 964          | 23           | 840 16/07/2793 03:49:23  | 23/01/2580 22:  |
| 6       | hNb4465q166cf315xU072ILq796uA86yQg686IG01IvZ      | 533          | 63           | 141 06/11/2494 07:04:48  | 28/05/2571 00:  |
| 7       | gQw888MB522vE915ySw395kP793U9318Y634Hu04kC        | 462          | 75           | 1795 27/02/2694 21:16:50 | 25/07/2224 07:  |
| 8       | yUn814kR289UA42kQF452vGH912pOL129Uw195jNB257N     | 973          | 2            | 957 06/01/2969 23:44:09  | 06/03/2707 04:  |
| 9       | PE691mGV184FT864K444PD572hMc457DQ780L187R         | 755          | 79           | 1205 26/09/2303 15:20:52 | 25/02/2625 03:  |
| 10      | FZA086uCC554oQH014oCc719Uq022d665dQy396hBo614V    | 390          | 79           | 1952 26/09/2019 03:58:45 | 28/01/2805 10:  |
| 11      | qZy727vGo437vWV390uJB160vY1120fC084kQ347wAY735vA  | 37           | 9            | 1807 06/03/2428 05:25:52 | 26/05/2513 18:  |
| 12      | kGU869hBU541hZ1216gVg797W6927FfX01QeV982bM494bR   | 784          | 7            | 266 20/06/2276 09:10:13  | 18/06/2914 08:  |
| 13      | vGq557zEo497AS116wWk259uJO352hGy181hU6393kAq562nS | 587          | 52           | 1238 04/09/2775 11:27:18 | 03/09/3018 03:  |
| 14      | cFW797aCV835oTc394gT4397Mh134Tt682eYf316Mg126vQ   | 618          | 25           | 1527 16/05/2119 03:02:00 | 12/06/2812 21:  |
| 15      | wf718aGe511jg889bBC090oO8563pIO357hJ741mIO15ZuK   | 662          | 59           | 662 19/06/2203 23:27:27  | 11/09/1921 01:  |
| 16      | dAB01GvV953yCw065uEJ1698h5459vY364nZ328hYL723vG   | 11           | 88           | 1985 19/10/2841 22:51:27 | 10/08/2428 01:  |
| 17      | vGH776hO26DhN587W784gEM559wLd090Yq727m3k144IG     | 447          | 89           | 2100 27/09/2379 03:16:42 | 15/07/2625 00:  |
| 18      | uLU698v6T71hKT469qQ756pNA350wF723oRq519vG549vK    | 971          | 44           | 2030 28/09/2788 16:14:18 | 07/11/2224 15:  |
| 19      | dUOS47aYr514uUj062vWD992vV025vYC344UB759kRD093K   | 133          | 35           | 1672 14/04/2381 14:46:01 | 21/10/2512 22:  |
| 20      | hVW621mKv318Rz22QY211PC789eYw428Uj646hRM789F      | 1217         | 6            | 1581 30/12/1954 11:54:23 | 13/06/3034 09:  |
| 21      | wY760Ph244bEm411aWq747hRe611uDC839Bm851Gz075IE    | 542          | 62           | 631 10/12/2624 02:09:01  | 04/09/2084 02:  |

1:1 Mailer@labdbwin 22 rows selected in 0.130 seconds (more...)

Enter filter text

Refrigerato

Stovel

Blood\_Donat

DeviceTest\_

Test\_Device

Test\_1

Test\_Tubel

## שאלות בחירה

1. שאלת שמחברת בין 3 טבלאות- מכשיר, סטט\_מכשיר וסטט. בשביל להוציא את פרטי הטסט. זה כולל תת שאלתא שסופרת את הסך הטסטים שנעשו עבור כל מכשיר.

```
SELECT
    d.DeviceId,
    d.DeviceType,
    t.TestID,
    t.TestDate,
    t.ResultT,
    (SELECT COUNT(*) FROM Test_Device1 td WHERE td.DeviceId = d.DeviceId) AS TotalTests
FROM
    Device1 d
JOIN
    Test_Device1 td ON d.DeviceId = td.DeviceId
JOIN
    Test_1 t ON td.TestID = t.TestID
ORDER BY
    d.DeviceId, t.TestDate;
```

בוחר את כל המכשירים עם הטסטים שלהם כולל את מספר הטסטים שנעשו על כל מכשיר--

```
SELECT
    d.DeviceId,
    d.DeviceType,
    t.TestID,
    t.TestDate,
    t.ResultT,
    (SELECT COUNT(*) FROM Test_Device1 td WHERE td.DeviceId = d.DeviceId) AS TotalTests
FROM
    Device1 d
JOIN
    Test_Device1 td ON d.DeviceId = td.DeviceId
JOIN
    Test_1 t ON td.TestID = t.TestID
ORDER BY
    d.DeviceId, t.TestDate;
```

|    | DEVICEID                | DEVICETYPE | TESTID                | TESTDATE   | RESULTT                | TOTALTESTS |
|----|-------------------------|------------|-----------------------|------------|------------------------|------------|
| 1  | 03waK1eR32N1aiv9OjMi    | 92         | 14UWhXp1P1jq7llyvniE  | 11/09/2006 | QF4vzIWCBdiaeDMcZpVR   | 2          |
| 2  | 03waK1eR32N1aiv9OjMi    | 92         | h9foIvpzzH1HNVkrT5p   | 28/06/2009 | suakQY7kSiUEvzyjeFTN   | 2          |
| 3  | 0o3lnEATBqky2vOiibhK    | 1          | 6dKTbnc00Ycrob45y5UK  | 21/10/2022 | kjqlqmrnfigwep8bTRVvk  | 1          |
| 4  | 18D0bmQF0pT08L9pAb6P    | 79         | nPK5Th5gZ0MvGfKox1    | 15/10/2008 | lclN34GutRZKv1bSnUGqD  | 3          |
| 5  | 18D0bmQF0pT08L9pAb6P    | 79         | flOzvnQH2JEfOZJ9Gnw   | 22/12/2008 | P6wo5yeh5pWenQXw1rpN   | 3          |
| 6  | 18D0bmQF0pT08L9pAb6P    | 79         | 6frcU4LsrnBgxfid0aU6  | 04/06/2021 | J9wZvgMkfd3dXlQ0r      | 3          |
| 7  | 1qMAFhgbPMb7p8U4ui      | 77         | mtHjZDjwycCkorpQnUn   | 19/04/2006 | EDFgl0BcqfK5jnj3pvh    | 2          |
| 8  | 1qMAFhgbPMb7p8U4ui      | 77         | lsmhwp5X9dY7Y35YR0pN  | 03/10/2017 | c8rpNLt8QD6Ma6YmcNCR   | 2          |
| 9  | 1suLrYJRcaaf0pdmM5g     | 70         | bxgyeq39EioMvvggEWbf  | 30/12/2008 | cwOpHdfpwb8dHvQDTxgUV  | 1          |
| 10 | 1FpJ0uEuB0omP8q3MBHw    | 72         | 8MuA9YvayfQCMMyE6Uj   | 01/01/2004 | B0V9fOpQic6GZPf7dfk    | 2          |
| 11 | 1K7C7uB6i17Dnm8BxcmABHw | 72         | 8W7G5EwO5A073w15A7P7C | 10/06/2013 | 47D77Tuu6K9d47M5CphmYm | 2          |

2. מחזירה את פרטי תרומת הדם יחד עם תוצאות הטסט ומחשבת את גיל תרומת הדם תוך שימוש בתאריך הנוכחי(כמה זמן עבר מאז שנתרם הדם).

```
SELECT
    bd.DonationID,
    bd.DonorID,
    bd.BloodType,
    bd.DonationDate,
    t.ResultT,
    TRUNC(SYSDATE - bd.DonationDate) AS AgeInDays
FROM
    Blood_Donation1 bd
JOIN
    Test_1 t ON bd.TestID = t.TestID
ORDER BY
    bd.DonationDate DESC;
```

-- מחזירה פרטים של תרומות דם ביחד עם תוצאות הטסט ומחשבת את גיל כל תרומות דם ביימים מהתאריך הנוכחי--  
-- מחזיר בסדר יורד של תאריכי טיס--

```
SELECT
    bd.DonationID,
    bd.DonorID,
    bd.BloodType,
    bd.DonationDate,
    t.ResultT,
    TRUNC(SYSDATE - bd.DonationDate) AS AgeInDays
FROM
    Blood_Donation1 bd
JOIN
    Test_1 t ON bd.TestID = t.TestID
ORDER BY
    bd.DonationDate DESC;
```

|    | BLOODTYPE  | DONATIONDATE        | RESULTT  | AGEINDAYS |
|----|--|---------------------|--|-----------|
| 2  | oOh147gBu285dxP071v8f558uGj3Q2fZ229qPv377kRf748pW  | 24/07/2064 03:22:01 | ObAw9aNNLJSgloQFvid                                | -379089   |
| 3  | lQP517uKd623aAh723zH112hZa896ePQ566Vl266TM859kY    | 14/09/2063 13:11:51 | 4OUA7mOj8k4kzZg7qJMb                               | -379575   |
| 4  | pCH788oYP298qYH477xG075yAg245cR021uOO289kH355aU    | 04/02/2063 02:23:16 | uxnUhhwCRv9kel3b8lM                                | -379953   |
| 5  | dk954ilI7xOY223Oj979wFO59k4DZ715zZl240yQC924tJ     | 16/07/2046 13:58:56 | v3v660gFN851uUX976gRd424mRz514yfN667dKM190bPH717u  | -373207   |
| 6  | yPv451lRj783qWAD43vTl894fRX628hZD942Cdm689vVR194pZ | 04/04/2043 06:27:21 | lDHAg6c7y60HfFDPkF6w                               | -372107   |
| 7  | fMR884u0169pWt669pMd067yAk381yEn237pIA153vOv524uJ  | 25/02/2017 20:31:20 | ynmu5W3yKultOPNHFqS                                | -362574   |
| 8  | fNS276Yv065yBC571LlJ310qYv643fLlA766FE860Ys765mW   | 01/06/2013 08:02:49 | 8xCSO25Mk178yNM8fZud                               | -361208   |
| 9  | uCU540Qz733oNz365IV039cCQ015WZ743Nm616gVv438fG     | 08/09/2010 23:20:15 | zLr968bOU871qV4529bGw647LAU770bOv827aQD895vSD546uB | -360212   |
| 10 | kFJ116Wx982dOw018lq421wBD325H4p656gTn477eOz687jM   | 10/08/2991 06:34:33 | qIB323mPb995hCW594oxY872bPL761Gd07FTc914cBI448yG   | -353243   |
| 11 | NW8446dU278hJt692aAH166eAD056oKc182Bd4660uPx216dM  | 10/02/2985 16:04:39 | DRAbITc8407pEffbxAX                                | -350872   |

3. מחזירה את פרטי מבחנת הטסט שכוללים את המיקום שלה ומחשבת את מספר הימים שותרו לה עד לתאריך התפוגה.

```
SELECT
    tt.TestTubeID,
    tt.LocationTestTube,
    tt.ExpiryDate,
    r.RefrigeratorId,
    r.LocationR,
    tt.ExpiryDate - SYSDATE AS DaysToExpiry
FROM
    Test_Tube1 tt
JOIN
    Refrigerator1 r ON tt.TestTubeID = r.TestTubeID
ORDER BY
    tt.ExpiryDate;
```

חזיר פרטי מבחנות הטסט, שכולל מיקום המיקום ומחשב את מספר הימים שנותרו עד לתאריך התפוגה--

```
SELECT
    tt.TestTubeID,
    tt.LocationTestTube,
    tt.ExpiryDate,
    r.RefrigeratorId,
    r.LocationR,
    tt.ExpiryDate - SYSDATE AS DaysToExpiry
FROM
    Test_Tube1 tt
JOIN
    Refrigerator1 r ON tt.TestTubeID = r.TestTubeID
ORDER BY
    tt.ExpiryDate;
```

| TESTTUBEID  | LOCATIONTESTTUBE     | EXPIRYDATE          | REFRIGERATORID                                     | LOCATIONR                                      |
|---|----------------------|---------------------|--|--|
| 1   xH145hDz297bWP928qHd497dDN693qEn715mCM473d6R494wf | USA                  | 30/04/1900 16:40:30 | a0O448Dv4311d4q340:8b299nH1692Oh777h4c768jU959pE   | LR1628Gv914xGL988HV085pHj158DA611tHY3778EP     |
| 2   BC099FR899CT332gLB97zeAm4265m136u8b027c5u913uR    | USA                  | 07/03/1922 22:44:24 | vPU1748fIN428pIQ213uSP415FZJ154nYg9791c278kYE714nQ | zBQ164vY629wP1264AI712w9W971qX0191jNp062pAZ    |
| 3   t7r554Gv216oKu740yVL189vYU987qBL390ZC0820mY570eR  | Italy                | 07/10/2016 10:59:13 | gWY265yMj670aTw240nQH799pIn667nTy491wRf687nPr968uN | xBL257JLT561mE:043yDL265Z2g042uZw597RvX983Uj59 |
| 4   t7r554Gv216oKu740yVL189vYU987qBL390ZC0820mY570eR  | Italy                | 07/10/2016 10:59:13 | wfU802cG:921aMG254eLUX5458m777Cb794Vz915kTb016eB   | fFU1548f782cVZ477kCF984yOP539wKQ935PfC465nPB   |
| 5   X3YnVa9z96uU95jnSzi                               | ChZefNOmorPX0nkDz8VA | 15/02/2024          | DDyLpBp93GLite8YKVF                                | iQDpVpYAvUuHIOF6p59m                           |

4. מחברת בין כמה טבלאות כדי לקבל פרטי מכשיר ותאריך שימוש אחרון. גם מחשבת את המספר הכללי של הטסטים שנעשו והזמן הממוצע שעבר ביניהם.

```
WITH TestIntervals AS (
    SELECT
        d.DeviceId,
        d.DeviceType,
        d.LastMaintenanceDate,
        t.TestID,
        t.TestDate,
        LAG(t.TestDate) OVER (PARTITION BY d.DeviceId ORDER BY t.TestDate) AS
PrevTestDate
    FROM
        Device1 d
    LEFT JOIN
        Test_Device1 td ON d.DeviceId = td.DeviceId
    LEFT JOIN
        Test_1 t ON td.TestID = t.TestID
)
SELECT
    DeviceId,
    DeviceType,
    LastMaintenanceDate,
    COUNT(TestID) AS NumberOfTests,
    AVG(TestDate - PrevTestDate) AS AvgTimeBetweenTests
FROM
    TestIntervals
GROUP BY
    DeviceId,
    DeviceType,
    LastMaintenanceDate
ORDER BY
    NumberOfTests DESC;
```

```
WITH TestIntervals AS (
    SELECT
        d.DeviceId,
        d.DeviceType,
        d.LastMaintenanceDate,
        t.TestID,
        t.TestDate,
        LAG(t.TestDate) OVER (PARTITION BY d.DeviceId ORDER BY t.TestDate) AS PrevTestDate
    FROM
        Device1 d
    LEFT JOIN
        Test_Device1 td ON d.DeviceId = td.DeviceId
    LEFT JOIN
        Test_1 t ON td.TestID = t.TestID
)
SELECT
    DeviceId,
    DeviceType,
    LastMaintenanceDate,
    COUNT(TestID) AS NumberOfTests,
    AVG(TestDate - PrevTestDate) AS AvgTimeBetweenTests
FROM
    TestIntervals
GROUP BY
    DeviceId,
    DeviceType,
    LastMaintenanceDate
```

|   | DEVICEID             | DEVICETYPE | LASTMAINTENANCEDATE | NUMBEROFTESTS | AVGTIMEBETWEENTESTS |
|---|----------------------|------------|---------------------|---------------|---------------------|
| 1 | DL0CQIDwQZOAHUFC86C  | ...        | 26 16/10/2017       | 5             | 82696.2433738426    |
| 2 | JfyxkJE9AWlpRA6FnVTh | ...        | 19 05/05/2005       | 5             | 95382.7063744213    |
| 3 | dk7GxRPa0LLsAi7kzNw4 | ...        | 49 05/11/2009       | 5             | 94483.7163975694    |
| 4 | VSGNUhHSO7effsHACTZ  | ...        | 96 27/03/2009       | 4             | 2602.6666666667     |

## שאלות עידכון

4. מעדכן מיקום למכשיר מסוג 83 שלא השתמשו בו למעלה משנה

5. UPDATE
6. Device1
7. SET
8. LocationD = 'Israel'
9. WHERE
10. DeviceType = 83
11. AND LastMaintenanceDate < ADD\_MONTHS(SYSDATE, -12);
- 12.

לפני:

```
SELECT* from Device1;
SELECT* from Test_Tube1;
SELECT* from Test_1;
SELECT* from Test_Device1;
SELECT* from DeviceTest_Tube1;
SELECT* from Blood_Donation1;
SELECT* from Stove1;
SELECT* from Refrigerator1;
```

| DEVIDEID | DEVICETYPE | SERIESNUM | VERSIOND | MANUFACTURER          | PURCHASEDATE | LASTMAINTENANCEDATE | LOCATIOND            |
|----------|------------|-----------|----------|-----------------------|--------------|---------------------|----------------------|
| 1        | 83         | 193       | 0        | 9E6uEKFOsINXZIRZVZm   | 17/06/2017   | 05/12/2012          | d79e4zWj/pnT6nu3uXBa |
| 2        | 31         | 760       | 9        | PvL2vp1kgRPSYTZaTE    | 16/09/2003   | 13/04/2010          | evqj5CvnUWT5sjpFOr8  |
| 3        | 46         | 436       | 2        | WfFvJLbwrlcYzZdpL5G   | 04/04/2010   | 25/09/2010          | mfFjXEDHUPMHZDlvApsT |
| 4        | 73         | 748       | 9        | sApvxdDVio2s9G4AMJ3C  | 12/03/2020   | 20/04/2000          | r964gYNL8H8739jYgrf  |
| 5        | 98         | 717       | 0        | SUOVQmny59Rv11gmYE9UT | 28/04/2015   | 17/09/2012          | 44liM4MUowLILqKcO3aF |
| 6        | 28         | 78        | 5        | kOoE53yAvrvNokKfoZ    | 21/07/2022   | 14/02/2022          | ZD0oIkOCpUln38j3TNS  |
| 7        | 97         | 566       | 0        | 4yz29FkoMl3yOo5LkO    | 02/01/2017   | 07/06/2021          | T2L8toOv80V7uBClBz4k |
| 8        | 72         | 454       | 9        | M6EHEDvYBTeE8NlJNGI   | 15/12/2003   | 16/06/2009          | nKlKSoHl1F5I79woA8b  |
| 9        | 81         | 945       | 5        | jqzW204LWuG9yZOKICIdS | 20/10/2014   | 17/06/2013          | wzhwN0gCbvLV5AashXz1 |
| 10       | 92         | 982       | 1        | jczdVIMgtOWi1bPOF0Pe  | 24/08/2005   | 22/09/2000          | Jr4ACbfnLDAmlea6aEx2 |
| 11       | 56         | 975       | 7        | wunWfIkIBHR7K0OZMg    | 15/10/2013   | 15/03/2001          | 6XU6fYgyvPINP27U6VOQ |
| 12       | 90         | 132       | 2        | D5Ui763SK6E65vK4H4I   | 04/12/2002   | 06/11/2000          | TY7m6wThJN67Z6x0MzqG |
| 13       | 16         | 424       | 5        | ZzGwKulUg7dQ7zrrdka   | 13/10/2015   | 28/04/2020          | KzH7jthfIUrmwz15LHu3 |

אחרי:

```
SELECT* from Device1;
SELECT* from Test_Tube1;
SELECT* from Test_1;
SELECT* from Test_Device1;
SELECT* from DeviceTest_Tube1;
SELECT* from Blood_Donation1;
SELECT* from Stove1;
SELECT* from Refrigerator1;
```

| DEVIDEID | DEVICETYPE | SERIESNUM | VERSIOND | MANUFACTURER          | PURCHASEDATE | LASTMAINTENANCEDATE | LOCATIOND            |
|----------|------------|-----------|----------|-----------------------|--------------|---------------------|----------------------|
| 1        | 83         | 193       | 0        | 9E6uEKFOsINXZIRZVZm   | 17/06/2017   | 05/12/2012          | Israel               |
| 2        | 31         | 760       | 9        | PvL2vp1kgRPSYTZaTE    | 16/09/2003   | 13/04/2010          | evqj5CvnUWT5sjpFOr8  |
| 3        | 46         | 436       | 2        | WfFvJLbwrlcYzZdpL5G   | 04/04/2010   | 25/09/2010          | mfFjXEDHUPMHZDlvApsT |
| 4        | 73         | 748       | 9        | sApvxdDVio2s9G4AMJ3C  | 12/03/2020   | 20/04/2000          | r964gYNL8H8739jYgrf  |
| 5        | 98         | 717       | 0        | SUOVQmny59Rv11gmYE9UT | 28/04/2015   | 17/09/2012          | 44liM4MUowLILqKcO3aF |

## 2. מעדכן שימוש אחרון ומיקום למכשירים שנעשה בהם שימוש בשלוש השנים האחרונות.

```

4. --update2
5.
6. מעדכן שימוש אחרון ומיקום למכשירים שנעשה בהם שימוש בשלוש השנים האחרונות--
7. commit
8. UPDATE
9. Device1
10. SET
11. LastMaintenanceDate = SYSDATE,
12. LocationD='USA'
13. WHERE
14. DeviceId IN (
15. SELECT td.DeviceId
16. FROM Test_Device1 td
17. JOIN Test_1 t ON td.TestID = t.TestID
18. WHERE t.TestDate > SYSDATE
19. );

```

לפני:

| DEVICEID | DEVICETYPE | SERIENUM | VERSIOND | MANUFACTURER                                      | PURCHASEDATE        | LASTMAINT           |
|----------|------------|----------|----------|---|---------------------|---------------------|
| 1        | 20         | 81       | 1        | vZv143e167aVU113dJU824mF5241vCw492mBH222gUL172hD  | 21/05/2549 05:04:42 | 30/07/2567 00:03:53 |
| 2        | 55         | 570      | 10       | vZT074mQj667yLw081Rn048vG6479ZU785yK5569Gf827aL   | 15/11/2386 20:06:17 | 29/02/2669 1        |
| 3        | 15         | 943      | 1        | hUs798vZ798vH677aOV714S1927Z3823Ao011bAw954dQ     | 18/11/1987 10:51:53 | 02/06/2709 0        |
| 4        | 6          | 765      | 10       | gZq198RZ41bOCT816pA432274679ZU115pK645fT21Z       | 17/01/2007 20:58:27 | 18/01/2618 1        |
| 5        | 60         | 927      | 6        | yA0211A4077ON893kM093wUY181wE6797mKc562cCF962wH   | 28/02/2045 05:18:05 | 14/05/2095 1        |
| 6        | 19         | 790      | 0        | Kxg830QfW951BK742jGQ342iNQT798uQIT46wY890bC492CZ  | 11/10/2216 09:23:28 | 24/01/2614 2        |
| 7        | 78         | 332      | 9        | oC540eQ672nRK397pYY460Xa218mOP520KxS472vV998nO    | 20/11/2443 00:01:35 | 15/03/2269 2        |
| 8        | 80         | 239      | 4        | jFC556eOK925dER039Nk433ST953Gn492oG4647cPE254v    | 21/05/2743 15:24:11 | 30/09/2502 1        |
| 9        | 66         | 8        | 0        | sSH424gUq132wUG496Uj348vA935pEx625Ym044nT665gl    | 07/09/2076 13:19:12 | 02/01/2924 1        |
| 10       | 31         | 321      | 9        | jCt065zAR128RH044Lj492gCj367dM5737uOU624ePM818gG  | 05/02/2416 17:18:57 | 14/07/2199 0        |
| 11       | 25         | 449      | 6        | xPK398Z0413xDc384Uj740dfg685hDR395bQJ757bSH971eE  | 24/02/2391 11:10:45 | 30/10/2682 0        |
| 12       | 51         | 995      | 1        | zKx887eUj438O1880kLw086uOy814Fxe691nUg956dPz962rB | 18/03/2620 14:25:35 | 27/06/2220 2        |
| 13       | 14         | 661      | 8        | oB1574IOV251kU587kV868LL180bOh258vOT021bSR950mF   | 21/07/2478 21:12:57 | 06/05/2709 1        |

אחרי:

| DEVICETYPE | SERIENUM | VERSIOND | MANUFACTURER | PURCHASEDATE                                       | LASTMAINTENANCEDATE | LOCATIOND |
|------------|----------|----------|--------------|--|---------------------|-----------|
| 1          | 20       | 81       | 1            | vZv143e167aVU113dJU824mF5241vCw492mBH222gUL172hD   | 21/05/2549 05:04:42 | USA       |
| 2          | 34       | 863      | 0            | rTg078aUk934kAp638Sh449wDZ455gHg590rAq535Fen614gL  | 04/09/2150 08:56:15 | USA       |
| 3          | 78       | 332      | 9            | sXs540eQf672nRK397pYY460Xa218mOP520KxS472vV998nO   | 20/11/2443 00:01:35 | USA       |
| 4          | 31       | 321      | 9            | jCt065zAR128RH044Lj492gCj367dM5737uOU624ePM818gG   | 05/02/2416 17:18:57 | USA       |
| 5          | 25       | 449      | 6            | xPK398Z0413xDc384Uj740dfg685hDR395bQJ757bSH971eE   | 24/02/2391 11:10:45 | USA       |
| 6          | 51       | 995      | 1            | zKx887eUj438O1880kLw086uOy814Fxe691nUg956dPz962rB  | 18/03/2620 14:25:35 | USA       |
| 7          | 14       | 661      | 8            | oB1574IOV251kU587kV868LL180bOh258vOT021bSR950mF    | 21/07/2478 21:12:57 | USA       |
| 8          | 21       | 29       | 2            | nWJ962dPm63DyFr164AX262oEZ615xXE186Ep517vWV753mN   | 10/07/2535 18:48:31 | USA       |
| 9          | 17       | 397      | 9            | rVv457ITTH692y1z117cTz855dV137wVY331mPK786cRX572d8 | 16/05/2147 23:37:06 | USA       |
| 10         | 41       | 656      | 5            | aVe877gPR619dFe031kY330oGq764pK622FQ141ITA355cXQ   | 12/04/2394 19:16:01 | USA       |
| 11         | 4        | 260      | 8            | pSu2341St645tOU250jKc696fBD682dSH778vOC127kV1887O  | 21/11/2549 12:41:19 | USA       |
| 12         | 5        | 877      | 6            | wRO572aMn430vY560dLr762Zzh364vVo827mMt6280Uj255uL  | 11/10/2989 08:07:35 | USA       |
| 13         | 85       | 144      | 8            | fRI152rDs167Jq460Vq652nLy639A3368YO649nD1983tG     | 02/10/2476 03:21:53 | USA       |



## שאלות מחיקה

4. מוחק תרומות דם שנעשו לפני יותר משנתיים. מוחק מכל הטבלאות הקשורות בזה.

```

5. BEGIN
6.     -- Delete related records in Refrigerator1
7.     DELETE FROM Refrigerator1
8.     WHERE DonationID IN (
9.         SELECT DonationID
10.        FROM Blood_Donation1
11.        WHERE DonationDate < ADD_MONTHS(SYSDATE, -24)
12.    );
13.
14.    -- Delete related records in Test_Device1
15.    DELETE FROM Test_Device1
16.    WHERE TestID IN (
17.        SELECT TestID
18.        FROM Test_1
19.        WHERE BloodDonationID IN (
20.            SELECT DonationID
21.            FROM Blood_Donation1
22.            WHERE DonationDate < ADD_MONTHS(SYSDATE, -24)
23.        )
24.    );
25.
26.    -- Delete related records in Blood_Donation1
27.    DELETE FROM Blood_Donation1
28.    WHERE DonationDate < ADD_MONTHS(SYSDATE, -24);
29.
30.    -- Delete records in Test_1
31.    DELETE FROM Test_1
32.    WHERE BloodDonationID IN (
33.        SELECT DonationID
34.        FROM Blood_Donation1
35.        WHERE DonationDate < ADD_MONTHS(SYSDATE, -24)
36.    );
37.

```

לפני:

```

SELECT* FROM Stove1;
SELECT* FROM Refrigerator1;

select * from Blood_Donation1
where DonationDate < ADD_MONTHS(SYSDATE, -24);

select * from Device1
where LastMaintenanceDate > SYSDATE
;
select * from Device1

```

|   | DONATIONID          | DONORID              | BLOODTYPE            | DONATIONDATE | VOLUMEDONATED | TESTRESULTID         | TESTID                  |
|---|---------------------|----------------------|----------------------|--------------|---------------|----------------------|-------------------------|
| 1 | Gy6Dz0mbrqChvSg7DqA | wwemDqUakyf7Phhw759z | KTZ6vFty3Cm4WkogV2q  | 14/02/2022   | 321           | qH10bcncpnjMdnXHVWz  | r1pTidWF7ujWgd26Gr8l    |
| 2 | kgIj278yESQP0NyZdpm | nPKgAZXTkuqY0sZelaa  | FXZHDnRShj7JUnSGap0g | 19/05/2022   | 328           | qXqg08eBYfoMySaVhQR4 | kV3waV0uujuuwOxs4jnYdQE |

אחרי:

```
SELECT * from Blood_Donation1
WHERE DonationDate < ADD_MONTHS(SYSDATE, -24);
SELECT * from Device1
WHERE VersionD =0;
```

| DONATIONID | DONORID | BLOODTYPE | DONATIONDATE | VOLUMEDONATED | TESTRESULTID | TESTID |
|------------|---------|-----------|--------------|---------------|--------------|--------|
|------------|---------|-----------|--------------|---------------|--------------|--------|

```
WHERE DeviceId IN (
SELECT DeviceId
FROM Device1
WHERE VersionD =0
);
```

```
-- Delete related records in Test_Device1
```

```
DELETE FROM Test_Device1
```

```
WHERE DeviceId IN (
SELECT DeviceId
FROM Device1
WHERE VersionD =0
);
```

```
-- Delete related records in DeviceTest_Tube1
```

```
DELETE FROM DeviceTest_Tube1
```

```
WHERE DeviceId IN (
SELECT DeviceId
FROM Device1
WHERE VersionD =0
);
```

```
-- Delete devices of type 83
```

```
DELETE FROM Device1
WHERE VersionD =0;
```

END;

לפני:

```
SELECT * from Device1
WHERE VersionD =0;
```

| DEVID | DEVICETYPE  | SERIESNUM | VERSIOND | MANUFACTURER   | PURCHASEDATE        | LASTMAINTEN     |
|-------|---|-----------|----------|--|---------------------|-----------------|
| 3     | GmrbtE8X9GMV5T3U16                                | 97        | 566      | 0 4yz28FKoM3yOo5LkO                                  | 02/01/2017          | 07/06/2021      |
| 4     | evFPXla550Tntc3wQA6                               | 81        | 356      | 0 gNLtYjAjrE1Dblma50s                                | 19/08/2007          | 11/09/2014      |
| 5     | grtD2XQAJgTZPtc6gRn                               | 49        | 975      | 0 iVe9yIsXh2Pg3MQOKN4N                               | 14/09/2004          | 17/02/2006      |
| 6     | dYCqiwEngdq1XkHDdw4                               | 52        | 252      | 0 O3RnDyYASoWcAJMmaFz                                | 24/12/2006          | 16/06/2024 15:4 |
| 7     | vssuOoN770YrOIlvsPVV                              | 78        | 346      | 0 PwWjuX8ydl3oEpGRNBrl                               | 21/04/2010          | 19/11/2007      |
| 8     | aRq239bWc744lPX032vYt722mOK220mG743jPz879bDi791rD | 34        | 863      | 0 rTg078aUk934kAp63jSh449wOZ455gHg590Aq535fEn614gL   | 04/09/2150 08:56:15 | 16/06/2024 15:4 |
| 9     | xMi989cSB761sHg465cAE340LD332yAw971gGR536aPz095oX | 19        | 790      | 0 kXq030oRw951tBK742cGQ342fNQ738uQ1744kVY890aBC492cZ | 11/10/2216 09:23:28 | 24/01/2614 23:2 |
| 10    | zUZ677LRK714RF180BX550EJ057hBU879Zy687eET829cT    | 66        | 8        | 0 sSH424gUq132wUg496lUj348vAA95pEz625Ym044eNT665gl   | 07/09/2076 13:19:12 | 02/01/2924 13:C |
| 11    | vFn171oTa222mWz244SY389cKD418pHr169qW0076ZV082xT  | 75        | 457      | 0 tMb596vWF817kQW323KL0976oJb750pCP498lIN124d1117vW  | 03/07/2786 10:56:38 | 16/06/2024 15:4 |
| 12    | zBy324vQv417INy039YD886cVz156dER439wDx334zSZ994bX | 67        | 815      | 0 kGX890vZ0971pJc557oRc283dL398bDE189cOo314nCV131rD  | 02/12/2013 08:02:41 | 09/03/2047 03:4 |
| 13    | uIZ657jPN064fP869wEL643nVr345vMH131vGI012rC5458pN | 6         | 974      | 0 kZM931mUP572qZB190nlq178mLa824mFT264wNy245dWj287vW | 15/10/2927 20:42:34 | 03/01/2657 11:1 |
| 14    | vOQ374kKp447c50261kCr298qQm140vO017aYb179bDG870XS | 25        | 598      | 0 uYJ128vSA335hHs482nZB658qNg627vY913bLV546vAr763wL  | 01/05/2481 12:05:36 | 30/12/2873 09:3 |

אחרי:

```
SELECT * from Device1
WHERE VersionD =0;
```


| DEVID | DEVICETYPE | SERIESNUM | VERSIOND | MANUFACTURER | PURCHASEDATE | LASTMAINTENANCE | LOCATIOND |
|-------|------------|-----------|----------|--------------|--------------|-----------------|-----------|
|-------|------------|-----------|----------|--------------|--------------|-----------------|-----------|

## שאלות עם פרמטרים

(כל הצילומי מסך הנצרכים נמצאים באותו צילום מסך...)

1. מחזיר פרטי מכשיר שנמצא במיקום מסוים וסוג המכשיר שלו נבחר מרשימת סוגי המכשירים האפשריים.

```
SELECT
    DeviceId,
    DeviceType,
    Manufacturer,
    LocationD,
    LastMaintenanceDate
FROM
    Device1
WHERE
    LocationD = '&location'
    AND DeviceType IN (&<name="devType" list="select Distinct DeviceType from
Device1 " multiselect="yes">)
```



The screenshot shows a database management tool interface. The top pane displays a SQL query that filters devices by location and device type. The bottom pane shows the results of this query as a table with five columns: DeviceID, DeviceType, Manufacturer, LocationD, and LastMaintenanceDate. Three rows of data are visible, each with a unique DeviceID, a DeviceType of 5, a Manufacturer, a LocationD of USA, and a LastMaintenanceDate.

|   | DeviceID  | DeviceType | Manufacturer                                       | LocationD | LastMaintenanceDate |
|---|---|------------|--|-----------|---------------------|
| 1 | ITi421nS538XU949VY540VG154FYF711aMt626hCN9234W    | 5          | wRO572aMn430bVY560dLr762zZh364V0827mM628oUj255uL   | USA       | 16/06/2024 15:46:52 |
| 2 | JWt780kRV740OU155eBm730aUJ381eUJ375oUF017bEm917jH | 85         | fRi152rDa167jng460Nvq652nLy639AU368YO649hDj983tG   | USA       | 12/11/1960 23:26:00 |
| 3 | zKp640oIG030vOu892d0144nCB555uMc939hSv713gXm078gY | 5          | bNM291uVX276o3b063qGT335HBE414uZV582vCK453rLR836vD | USA       | 07/07/2363 02:20:58 |

## 2. מחזיר פרטי מבחנת סטט שעומדת בקיבול שניתנה ובתאריך תפוגה שהוקש.

```
SELECT
    TestTubeID,
    Material,
    VolumeCapacity,
    ExpiryDate,
    LocationTestTube
FROM
    Test_Tube1
WHERE
    VolumeCapacity > &<name="volumeCapacity">
    AND ExpiryDate < &<name="Expdate" type="date" required=true>
ORDER BY
    VolumeCapacity DESC;
```

The screenshot displays a database application interface. The top panel shows a SQL query editor with the following query:

```
SELECT
    TestTubeID,
    Material,
    VolumeCapacity,
    ExpiryDate,
    LocationTestTube
FROM
    Test_Tube1
WHERE
    VolumeCapacity > &<name="volumeCapacity">
    AND ExpiryDate < &<name="Expdate" type="date" required=true>
ORDER BY
    VolumeCapacity DESC;
```

A 'Variables' dialog box is open, showing the following values:


| Name           | Value      |
|----------------|------------|
| volumeCapacity | 120        |
| Expdate        | 01/01/2030 |

The bottom panel shows the results of the query in a table with 5 columns: TESTTUBEID, MATERIAL, VOLUMECAPACITY, EXPIRYDATE, and LOCATIONTESTTUBE. The table contains 13 rows of data.

| TESTTUBEID | MATERIAL             | VOLUMECAPACITY | EXPIRYDATE     | LOCATIONTESTTUBE     |
|------------|----------------------|----------------|----------------|----------------------|
| 1          | BPx7Q5ajOgDyaymXyyhe | 56             | 499 10/01/2025 | zu0FHfeGPqroleRNF1J8 |
| 2          | u2SDQnQDWAfUe6FutlDs | 56             | 498 11/02/2029 | FDB2wfFje1o9dnnMn33  |
| 3          | FsiY90/VKToZtgEnGCao | 40             | 493 26/06/2014 | eX8q749ZIHhtMMID:Lp7 |
| 4          | Vfszel8FequQllp0uLKg | 91             | 493 25/04/2009 | WVP2YygWnSjNGUAArObf |
| 5          | Vhh703ztP1tZ7GoZbWx  | 37             | 491 12/10/2024 | HkoPkerdANrHxUBECuA  |
| 6          | IgEzW2nTs8uUsYpC2tDe | 82             | 491 21/09/2029 | ciRDbfMZWW1Jee8yYEA  |
| 7          | qPLMLhpuSTGG8k4w6CW  | 92             | 487 25/03/2025 | VY21eCtpXsrMjL43xG9Z |
| 8          | fw8vScVybteAl6yZY2D9 | 56             | 486 08/10/2021 | yUKIANOBbuCKJBTlNeoB |
| 9          | aiV9KjwbE6F3pmLV91pR | 96             | 484 21/06/2029 | c8YuhV2mkILuZGY5PrQ  |
| 10         | JISjE5nSCf73HoP09mCG | 58             | 483 02/09/2015 | vsVuMpoXo3Qt15wNUj9  |
| 11         | mO9gBCFifvPNlWAh64De | 99             | 482 26/05/2010 | aO6hvaHacWlyZ6N7wJ3d |
| 12         | Vdq9FHHtOalQhwZQfa1S | 38             | 479 20/03/2025 | gi3OqftRmMqVGDruhCuR |
| 13         | nZCfK3UFFO746eW4FyN  | 4              | 478 26/05/2028 | PEcNf8UunnO8BfoLL2Gg |


3. מחזיר פרטי טסט עבור תוצאה שהתקבלה ושנמצא בין 2 תאריכים שהוקשו. התאריך סיום הינו אופציונלי.

```
SELECT
    TestID,
    MainTestType,
    SubTestType,
    TestDate,
    ResultT,
    BloodDonationID
FROM
    Test_1
WHERE
    ResultT = &<name="result" type="string">
    AND TestDate BETWEEN &<name="startDate" type="date" required=true> AND
    &<name="endDate" type="date" default="01/01/2030">
ORDER BY
    TestDate DESC;
```



4. מחזיר פרטי תרומת דם כאשר הID שלה נבחר מבין רשימת ID's וסוג הדם מוקש.

```
SELECT
    DonationID,
    DonorID,
    BloodType,
    DonationDate,
    VolumeDonated,
    TestResultID
FROM
    Blood_Donation1
WHERE
    DonorID IN &<name="donorIds"
list="'4vBdG7waNm9gJnMXKQLx','1CB675rNb061bRN689iPf743jMj288lMQ733nBi871sFR059v
A','1QC458wWb454lIk918eKK867mXJ138cTQ765yBq713qKp529eA'">
    AND BloodType = &<name="bloodType" type="string">
ORDER BY
    DonationDate DESC;
```




## אילוצים

1. הוספת אילוץ check ל-TestTube שבדוק שמוכנס רק ערך חיובי ל-VolumeCapacity.

```
ALTER TABLE Test_Tube1
ADD CONSTRAINT chk_VolumeCapacity CHECK (VolumeCapacity > 0);

INSERT INTO Test_Tube1 (TestTubeID, Material, VolumeCapacity, PurchaseDate, ExpiryDate)
VALUES ('TT001', 1, -10, TO_DATE('2024-01-01', 'YYYY-MM-DD'), TO_DATE('2025-01-01', 'YYYY-MM-DD'));
```



Error

ORA-02290: check constraint (OROKACH.CHK\_VOLUMECAPACITY) violated


OK Cancel Help

2. הוספת אילוץ ל-BLOOD DONATION שעבורו השדה DONORID יהיה ייחודי.

```
ALTER TABLE Blood_Donation1
ADD CONSTRAINT unique_DonorID UNIQUE (DonorID);

INSERT INTO Blood_Donation1 (DonationID, DonorID, BloodType, DonationDate, VolumeDonated, TestResultID, TestID)
VALUES ('BD001', 'DON001', 'A+', TO_DATE('2023-01-01', 'YYYY-MM-DD'), 500, 'TRO01', 'T001');

-- Attempt to insert a duplicate DonorID
INSERT INTO Blood_Donation1 (DonationID, DonorID, BloodType, DonationDate, VolumeDonated, TestResultID, TestID)
VALUES ('BD002', 'DON001', 'B+', TO_DATE('2023-02-01', 'YYYY-MM-DD'), 450, 'TRO02', 'T002');
```



Error

ORA-02291: integrity constraint (OROKACH.SYS\_C00708277) violated - parent key not found


OK Cancel Help

3. הוספת ערך ברירת מחדל ל-Device1 עבור Manufacturer.

```
ALTER TABLE Device1
MODIFY Manufacturer DEFAULT 'Unknown';

INSERT INTO Device1 (DeviceID, DeviceType, SeriesNum, VersionD, PurchaseDate, LastMaintenanceDate, LocationD)
VALUES ('D0001', 1, 100, 1, TO_DATE('2023-01-01', 'YYYY-MM-DD'), TO_DATE('2023-06-01', 'YYYY-MM-DD'), 'Lab1');
```

```
select * from Device1
where Manufacturer = 'Unknown';
```



|   | DEVICEID | DEVICETYPE | SERIESNUM | VERSIOND | MANUFACTURER | PURCHASEDATE | LASTMAINTENANCEDATE | LOCATIOND  |
|---|----------|------------|-----------|----------|--------------|--------------|---------------------|------------|
| 1 | D0001    | ...        | 1         | 100      | 1            | Unknown      | ...                 | 01/01/2023 |
|   |          |            |           |          |              |              |                     | 01/06/2023 |
|   |          |            |           |          |              |              |                     | Lab1       |

## Procedures

1. הוספת מכשיר חדש לבסיס הנתונים בכל מקרה שחסר נתון מסוים, זורקת חריגה שאומרת איזה נתון חסר.

```
) CREATE OR REPLACE FUNCTION mostDonations
    ,FromDate IN DATE
    ,ToDate IN DATE
    RETURN VARCHAR2 (
        IS
        ;(100)FunctionResult VARCHAR2
        ;MostDonationCount NUMBER
        ;(50)MostDonorID VARCHAR2
        BEGIN
        BEGIN
        SELECT DonorID, COUNT(DonationID) AS DonationCount
            INTO MostDonorID, MostDonationCount
            FROM Blood_Donation1
        WHERE DonationDate BETWEEN FromDate AND ToDate
            GROUP BY DonorID
            ORDER BY DonationCount DESC
        FETCH FIRST 1 ROWS ONLY; -- Use the correct Oracle syntax for fetching the
                                   first row

        ;FunctionResult := MostDonorID || ', ' || MostDonationCount
        EXCEPTION
        WHEN NO_DATA_FOUND THEN
        ;FunctionResult := 'No donations found in the given date range
        WHEN OTHERS THEN
        ;FunctionResult := 'An error occurred: ' || SQLERRM
        ;END
        ;RETURN FunctionResult
        ;END mostDonations
```



```

) CREATE OR REPLACE FUNCTION validateMaintenance
    AllowedPeriod IN NUMBER
    RETURN VARCHAR2 (
        IS
        CURSOR device_cursor IS
        SELECT DeviceId, LastMaintenanceDate
        ;FROM Device1

device_rec device_cursor%ROWTYPE;--for calculating dates
;e_OverdueMaintenance EXCEPTION
;(200)v_error_message VARCHAR2
;overdue_days NUMBER

BEGIN
;OPEN device_cursor

LOOP
;FETCH device_cursor INTO device_rec
;EXIT WHEN device_cursor%NOTFOUND

overdue_days := TRUNC(SYSDATE) -
;TRUNC(device_rec.LastMaintenanceDate)

Check if the device's last maintenance date exceeds the allowed period --
AllowedPeriod THEN < IF overdue_days
v_error_message := 'Device ' || device_rec.DeviceId || ' is overdue for
;.maintenance by ' || (overdue_days - AllowedPeriod) || ' days
;RAISE e_OverdueMaintenance
;END IF
;END LOOP

;CLOSE device_cursor

;RETURN 'All devices are within the allowed maintenance period

EXCEPTION
WHEN e_OverdueMaintenance THEN
;RETURN v_error_message
WHEN OTHERS THEN
;RETURN 'An unexpected error occurred: ' || SQLERRM
;END validateMaintenance

```

```

)CREATE OR REPLACE PROCEDURE InsertNewDevice
    ,IN VARCHAR2      p_DeviceId
    ,IN NUMBER        p_DeviceType
    ,IN NUMBER        p_SeriesNum
    ,IN NUMBER        p_VersionD
    ,IN VARCHAR2      p_Manufacturer
    ,IN DATE          p_PurchaseDate
    ,p_LastMaintenanceDate IN DATE
    IN VARCHAR2      p_LocationD
    IS (
    Define custom exception for missing values --
        ;e_MissingValue EXCEPTION
        ;(100)v_MissingField VARCHAR2
        BEGIN
            Check for missing values --
            IF p_DeviceId IS NULL THEN
                ;'v_MissingField := 'DeviceId
                ;RAISE e_MissingValue
            ELSIF p_DeviceType IS NULL THEN
                ;'v_MissingField := 'DeviceType
                ;RAISE e_MissingValue
            ELSIF p_SeriesNum IS NULL THEN
                ;'v_MissingField := 'SeriesNum
                ;RAISE e_MissingValue
            ELSIF p_VersionD IS NULL THEN
                ;'v_MissingField := 'VersionD
                ;RAISE e_MissingValue
            ELSIF p_Manufacturer IS NULL THEN
                ;'v_MissingField := 'Manufacturer
                ;RAISE e_MissingValue
            ELSIF p_PurchaseDate IS NULL THEN
                ;'v_MissingField := 'PurchaseDate
                ;RAISE e_MissingValue
            ELSIF p_LastMaintenanceDate IS NULL THEN
                ;'v_MissingField := 'LastMaintenanceDate
                ;RAISE e_MissingValue
            ELSIF p_LocationD IS NULL THEN
                ;'v_MissingField := 'LocationD
                ;RAISE e_MissingValue
            ;END IF

            Insert the new device record --
            INSERT INTO Device1 (DeviceId, DeviceType, SeriesNum, VersionD,
                Manufacturer, PurchaseDate, LastMaintenanceDate, LocationD)
                VALUES (p_DeviceId, p_DeviceType, p_SeriesNum, p_VersionD,
;p_Manufacturer, p_PurchaseDate, p_LastMaintenanceDate, p_LocationD)


;COMMIT

```


```

EXCEPTION
    WHEN e_MissingValue THEN
;DBMS_OUTPUT.PUT_LINE('Error: Missing value for ' || v_MissingField)
;ROLLBACK
    WHEN OTHERS THEN
;DBMS_OUTPUT.PUT_LINE('An unexpected error occurred: ' || SQLERRM)
;ROLLBACK
;END InsertNewDevice


```



לאחר הרצה:



חריגה:



2. עבור שני תז של מכשירים מחזיר רשימה לפי תאריכים עבור הימים שהיה להם את אותו מספר טסטים.

```

) CREATE OR REPLACE PROCEDURE FindSameTestDate
    ,p_DeviceId1 IN VARCHAR2
    ,p_DeviceId2 IN VARCHAR2

```

```

IS (
    Variables to store test counts --
    ;v_Device1_TestCount INT := 0
    ;v_Device2_TestCount INT := 0
    ;v_TestDate DATE

    Cursor to fetch test counts for each device on each date --
    CURSOR test_cursor IS
    ,SELECT t.TestDate
COUNT(CASE WHEN td.DeviceId = p_DeviceId1 THEN 1 END) AS
    ,Device1_TestCount
COUNT(CASE WHEN td.DeviceId = p_DeviceId2 THEN 1 END) AS
    ,Device2_TestCount
FROM Test_Device1 td
JOIN Test_1 t ON td.TestID = t.TestID
GROUP BY t.TestDate
HAVING COUNT(CASE WHEN td.DeviceId = p_DeviceId1 THEN 1 END) =
    ;COUNT(CASE WHEN td.DeviceId = p_DeviceId2 THEN 1 END)

    BEGIN
    Open the cursor --
    ;OPEN test_cursor

    Loop through the cursor results --
    LOOP
    FETCH test_cursor INTO v_TestDate, v_Device1_TestCount,
    ;v_Device2_TestCount
    ;EXIT WHEN test_cursor%NOTFOUND

    Compare and print results when counts are equal --
    DBMS_OUTPUT.PUT_LINE('On ' || TO_CHAR(v_TestDate, 'YYYY-MM-DD') || ', '
    || p_DeviceId1 || ' had ' || v_Device1_TestCount || ' tests and ' || p_DeviceId2 || ' had '
    ;|| v_Device2_TestCount || ' tests.')
    ;END LOOP

    Close the cursor --
    ;CLOSE test_cursor

    EXCEPTION
    WHEN NO_DATA_FOUND THEN
    DBMS_OUTPUT.PUT_LINE('No common test dates found for ' || p_DeviceId1 || '
    ;and ' || p_DeviceId2)
    WHEN OTHERS THEN
    ;DBMS_OUTPUT.PUT_LINE('An error occurred: ' || SQLERRM)
    ;END FindSameTestDate

```

077

| <input type="checkbox"/>            | Variable    | Type   | Value  |
|-------------------------------------|-------------|--------|--|
| <input checked="" type="checkbox"/> | p_deviceid1 | String | jTV718uPr846wIB933gWx760fAu549tVI320kFQ477sYI616eB |
| <input checked="" type="checkbox"/> | p_deviceid2 | String | fEs595bQe471kGH987aPw743nYD230zMt625pUs867qGv166oK |
| <input checked="" type="checkbox"/> | *           |        |  |

## פלט ההרצה:

```

On 2021-03-20, jTV718uPr846wIB933gWx760fAu549cVI320kFQ477sY1616eB had 0 tests and fEa595bQe471kG8987aPw743nYD230zM1625pUs867qOv166oK had 0 tests.
On 2021-04-26, jTV718uPr846wIB933gWx760fAu549cVI320kFQ477sY1616eB had 0 tests and fEa595bQe471kG8987aPw743nYD230zM1625pUs867qOv166oK had 0 tests.
On 2005-06-03, jTV718uPr846wIB933gWx760fAu549cVI320kFQ477sY1616eB had 0 tests and fEa595bQe471kG8987aPw743nYD230zM1625pUs867qOv166oK had 0 tests.
On 2008-04-25, jTV718uPr846wIB933gWx760fAu549cVI320kFQ477sY1616eB had 0 tests and fEa595bQe471kG8987aPw743nYD230zM1625pUs867qOv166oK had 0 tests.
On 2009-06-28, jTV718uPr846wIB933gWx760fAu549cVI320kFQ477sY1616eB had 0 tests and fEa595bQe471kG8987aPw743nYD230zM1625pUs867qOv166oK had 0 tests.
On 2022-07-10, jTV718uPr846wIB933gWx760fAu549cVI320kFQ477sY1616eB had 0 tests and fEa595bQe471kG8987aPw743nYD230zM1625pUs867qOv166oK had 0 tests.
On 2010-12-01, jTV718uPr846wIB933gWx760fAu549cVI320kFQ477sY1616eB had 0 tests and fEa595bQe471kG8987aPw743nYD230zM1625pUs867qOv166oK had 0 tests.
On 2013-03-19, jTV718uPr846wIB933gWx760fAu549cVI320kFQ477sY1616eB had 0 tests and fEa595bQe471kG8987aPw743nYD230zM1625pUs867qOv166oK had 0 tests.
On 2008-04-16, jTV718uPr846wIB933gWx760fAu549cVI320kFQ477sY1616eB had 0 tests and fEa595bQe471kG8987aPw743nYD230zM1625pUs867qOv166oK had 0 tests.
On 2020-11-14, jTV718uPr846wIB933gWx760fAu549cVI320kFQ477sY1616eB had 0 tests and fEa595bQe471kG8987aPw743nYD230zM1625pUs867qOv166oK had 0 tests.
On 2003-03-10, jTV718uPr846wIB933gWx760fAu549cVI320kFQ477sY1616eB had 0 tests and fEa595bQe471kG8987aPw743nYD230zM1625pUs867qOv166oK had 0 tests.
On 2002-09-04, jTV718uPr846wIB933gWx760fAu549cVI320kFQ477sY1616eB had 0 tests and fEa595bQe471kG8987aPw743nYD230zM1625pUs867qOv166oK had 0 tests.
On 2003-04-21, jTV718uPr846wIB933gWx760fAu549cVI320kFQ477sY1616eB had 0 tests and fEa595bQe471kG8987aPw743nYD230zM1625pUs867qOv166oK had 0 tests.
On 2006-12-07, jTV718uPr846wIB933gWx760fAu549cVI320kFQ477sY1616eB had 0 tests and fEa595bQe471kG8987aPw743nYD230zM1625pUs867qOv166oK had 0 tests.
On 3022-04-25, jTV718uPr846wIB933gWx760fAu549cVI320kFQ477sY1616eB had 0 tests and fEa595bQe471kG8987aPw743nYD230zM1625pUs867qOv166oK had 0 tests.
On 2017-07-21, jTV718uPr846wIB933gWx760fAu549cVI320kFQ477sY1616eB had 0 tests and fEa595bQe471kG8987aPw743nYD230zM1625pUs867qOv166oK had 0 tests.
On 2012-06-01, jTV718uPr846wIB933gWx760fAu549cVI320kFQ477sY1616eB had 0 tests and fEa595bQe471kG8987aPw743nYD230zM1625pUs867qOv166oK had 0 tests.

```

## Functions

1. פונקציה שמחזירה תז וכמות תרומות של אדם שתרם הכי הרבה תרומות דם בין שני תאריכים שניתנו כפרמטרים.

### הפונקציה:

```

) CREATE OR REPLACE FUNCTION mostDonations
    ,FromDate IN DATE
    ,ToDate IN DATE
    RETURN VARCHAR2 (
        IS
        ;(100)FunctionResult VARCHAR2
        ;MostDonationCount NUMBER
        ;(50)MostDonorID VARCHAR2
        BEGIN
        BEGIN
        SELECT DonorID, COUNT(DonationID) AS DonationCount
            INTO MostDonorID, MostDonationCount
            FROM Blood_Donation1
        WHERE DonationDate BETWEEN FromDate AND ToDate
            GROUP BY DonorID
            ORDER BY DonationCount DESC
        FETCH FIRST 1 ROWS ONLY; -- Use the correct Oracle syntax for fetching the
                                   first row

        ;FunctionResult := MostDonorID || ', ' || MostDonationCount
        EXCEPTION
        WHEN NO_DATA_FOUND THEN
        ;FunctionResult := 'No donations found in the given date range
        WHEN OTHERS THEN
        ;FunctionResult := 'An error occurred: ' || SQLERRM
        ;END
        ;RETURN FunctionResult
        ;END mostDonations
    
```

### התוצאה:

| Variable | Type   | Value   |
|----------|--------|---|
| result   | String | gTM888hFN991mRB268jBB036sEQ011bKu550gWZ591zEF796tF, 1 |
| fromdate | Date   | 02/05/2060  |
| todate   | Date   | 02/05/2090  |

2. עבור מספר ימים שמתקבל, מחזירה את המכשיר הראשון שמאז תאריך התחזוקה האחרון שלו עברו יותר ימים.

### הפונקציה:

```

) CREATE OR REPLACE FUNCTION validateMaintenance
    AllowedPeriod IN NUMBER
    RETURN VARCHAR2 (
        IS
        CURSOR device_cursor IS
        SELECT DeviceId, LastMaintenanceDate
        ;FROM Device1

device_rec device_cursor%ROWTYPE;--for calculating dates
;e_OverdueMaintenance EXCEPTION
;(200)v_error_message VARCHAR2
;overdue_days NUMBER

BEGIN
;OPEN device_cursor

LOOP
;FETCH device_cursor INTO device_rec
;EXIT WHEN device_cursor%NOTFOUND

overdue_days := TRUNC(SYSDATE) -
;TRUNC(device_rec.LastMaintenanceDate)

Check if the device's last maintenance date exceeds the allowed period --
AllowedPeriod THEN < IF overdue_days
v_error_message := 'Device ' || device_rec.DeviceId || ' is overdue for
;'.maintenance by ' || (overdue_days - AllowedPeriod) || ' days
;RAISE e_OverdueMaintenance
;END IF
;END LOOP

;CLOSE device_cursor

;RETURN 'All devices are within the allowed maintenance period

EXCEPTION
WHEN e_OverdueMaintenance THEN
;RETURN v_error_message
WHEN OTHERS THEN
;RETURN 'An unexpected error occurred: ' || SQLERRM
;END validateMaintenance

```

### התוצאה:

| Variable      | Type   | Value   |
|---------------|--------|---|
| result        | String | Device IT421nSk538sXu948vYs540vSI154FYF711aMt626hCN923tW is overdue for maintenance by 18 days. |
| allowedperiod | Float  | 10  |
|               |        |   |

## Mains

.1

התכנית:

```
CREATE OR REPLACE PROCEDURE MainProcedure
IS
    Parameters for FindSameTestDate --
    ;'v_DeviceId1 VARCHAR2(50) := 'D001
    ;'v_DeviceId2 VARCHAR2(50) := 'D002

    Parameters for validateMaintenance --
    v_AllowedPeriod NUMBER := 180; -- Example allowed period in days for
                                     maintenance
    ;(200)v_ValidationResult VARCHAR2

    BEGIN
        Call the FindSameTestDate procedure --
        BEGIN
            )FindSameTestDate
            ,v_DeviceId1 <= p_DeviceId1
            v_DeviceId2 <= p_DeviceId2
            ;(
            EXCEPTION
            WHEN OTHERS THEN
                DBMS_OUTPUT.PUT_LINE('An error occurred while calling FindSameTestDate:
                                     ;' || SQLERRM)
                ;END

        Call the validateMaintenance function --
        BEGIN
            ;v_ValidationResult := validateMaintenance(v_AllowedPeriod)
            DBMS_OUTPUT.PUT_LINE('Result of validateMaintenance: ' ||
                                ;v_ValidationResult)
            EXCEPTION
            WHEN OTHERS THEN
                DBMS_OUTPUT.PUT_LINE('An error occurred while calling
                                     ;validateMaintenance: ' || SQLERRM)
                ;END

        ;END MainProcedure
```



## התכנית:

```

CREATE OR REPLACE PROCEDURE MainProcedure
IS
    v_FromDate DATE := TO_DATE('2024-01-01', 'YYYY-MM-DD'); -- Example start
                                                                date
    Example end date -- ;v_ToDate DATE := TO_DATE('2024-12-31', 'YYYY-MM-DD')
                                                                ;(100)v_Result VARCHAR2

    Parameters for InsertNewDevice procedure --
        ;v_DeviceId VARCHAR2(50) := 'D001
        ;v_DeviceType NUMBER := 1
        ;v_SeriesNum NUMBER := 123
        ;v_VersionD NUMBER := 1
        ;v_Manufacturer VARCHAR2(50) := 'ManufacturerA
        ;v_PurchaseDate DATE := TO_DATE('2024-07-01', 'YYYY-MM-DD')
        ;v_LastMaintenanceDate DATE := TO_DATE('2024-07-01', 'YYYY-MM-DD')
        ;v_LocationD VARCHAR2(50) := 'LocationA
    BEGIN
        Call the mostDonations function --
        BEGIN
            ;v_Result := mostDonations(v_FromDate, v_ToDate)
            ;DBMS_OUTPUT.PUT_LINE('Result of mostDonations: ' || v_Result)
        EXCEPTION
            WHEN OTHERS THEN
                DBMS_OUTPUT.PUT_LINE('An error occurred while calling mostDonations: ' ||
                                                                ;SQLERRM)
                ;END

        Call the InsertNewDevice procedure --
        BEGIN
            )InsertNewDevice
            ,v_DeviceId <=      p_DeviceId
            ,v_DeviceType <=    p_DeviceType
            ,v_SeriesNum <=     p_SeriesNum
            ,v_VersionD <=      p_VersionD
            ,v_Manufacturer <=  p_Manufacturer
            ,v_PurchaseDate <=   p_PurchaseDate
            ,v_LastMaintenanceDate <= p_LastMaintenanceDate
            v_LocationD <=      p_LocationD
            ;(
            EXCEPTION
            WHEN OTHERS THEN
                DBMS_OUTPUT.PUT_LINE('An error occurred while calling InsertNewDevice: '
                                                                ;|| SQLERRM)
                ;END
        END MainProcedure
        ;

```

בסיס הנתונים לאחר הרצה:

